



CGNL1 rabbit pAb

Catalog No	YP-Ab-09316
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CGNL1 JACOP KIAA1749
Immunogen	Synthesized peptide derived from human CGNL1 AA range: 394-444
Specificity	This antibody detects endogenous levels of human CGNL1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CGNL1;CGNL1 JACOP KIAA1749;Cingulin-like protein 1 (Junction-associated coiled-coil protein) (Paracingulin)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell junction, tight junction . Localizes to the apical junction complex composed of tight and adherens junctions (PubMed:22891260). In the liver and kidney, it is also found along non-junctional actin filament bundles in addition to the apical junction (By similarity). .
Tissue Specificity	Smooth muscle, spleen, testis, fetal brain, amygdala, corpus callosum, cerebellum, thalamus and subthalamic nucleus of adult brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images