



CX04A rabbit pAb

Catalog No	YP-Ab-09317
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CXorf40A CXorf40 EOLA1
Immunogen	Synthesized peptide derived from human CX04A AA range: 15-65
Specificity	This antibody detects endogenous levels of human CX04A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CX04A;CXorf40A CXorf40 EOLA1;Protein CXorf40A (Endothelial-overexpressed lipopolysaccharide-associated factor 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed primarily in heart, skeletal muscle, kidney, liver and placenta. Relatively high level of expression in spleen, colon and small intestine. Almost no expression in brain, thymus, lung and peripheral blood leukocytes. Expressed in epithelial cells (at protein level) (PubMed:31007603).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images