



HS3S1 rabbit pAb

Catalog No	YP-Ab-09318
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HS3ST1 3OST 3OST1
Immunogen	Synthesized peptide derived from human HS3S1 AA range: 201-251
Specificity	This antibody detects endogenous levels of human HS3S1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	3-OST-1;3OST;3OST1;EC 2.8.2.23;h3 OST 1;h3-OST-1;Heparan sulfate (glucosamine) 3 O sulfotransferase 1;Heparan sulfate 3 O sulfotransferase 1; Heparan sulfate 3-O-sulfotransferase 1;Heparan sulfate D glucosaminyl 3 O sulfotransferase 1;Heparan sulfate D-glucosaminyl 3-O-sulfotransferase 1; Heparan sulfate glucosamine 3-O-sulfotransferase 1;Heparin glucosamine 3 O sulfotransferase;HS3S1_HUMAN;HS3ST 1;Hs3st1;HS3S1;HS3ST1 3OST 3OST1;Heparan sulfate glucosamine 3-O-sulfotransferase 1 (EC 2.8.2.23) (Heparan sulfate D-glucosaminyl 3-O-sulfotransferase 1) (3-OST-1) (Heparan sulfate 3-O-sulfotransferase 1) (h3-OST-1)
Observed Band	
Calculated Molecular Weight	35-5kD; 36kD
Cell Pathway	Golgi apparatus lumen .
Tissue Specificity	Highly expressed in the brain and kidney and weakly expressed in the heart, lung and placenta.
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

