



CCD34 rabbit pAb

Catalog No	YP-Ab-09325
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CCDC34 RAMA3 L15
Immunogen	Synthesized peptide derived from human CCD34 AA range: 249-299
Specificity	This antibody detects endogenous levels of human CCD34
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCD34_HUMAN;Ccdc34;Coiled-coil domain-containing protein 34;L15;NY-REN-41;RAMA3;Renal carcinoma antigen NY-REN-41;CCD34;CCDC34 RAMA3 L15;Coiled-coil domain-containing protein 34 (Renal carcinoma antigen NY-REN-41)
Observed Band	
Calculated Molecular Weight	43kDa, 35 kDa; 43kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

