



TYB4 rabbit pAb

Catalog No	YP-Ab-09327
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMSB4X TB4X THYB4 TMSB4
Immunogen	Synthesized peptide derived from human TYB4 AA range: 1-24
Specificity	This antibody detects endogenous levels of human TYB4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Fx;Hematopoietic system regulatory peptide;Prothymosin beta 4;PTMB 4; PTMB4;Seraspenide;T beta 4;T beta-4;TB4X;THYB 4;Thyb4;Thymosin beta 4; Thymosin beta 4 X chromosome;Thymosin beta 4 X linked;TMSB 4;TMSB4; Tmsb4x;TYB4 HUMAN;TYB4;TMSB4X TB4X THYB4 TMSB4;Thymosin beta-4 (T beta-4) (Fx) [Cleaved into: Hematopoietic system regulatory peptide (Seraspenide)]
Observed Band	
Calculated Molecular Weight	5 kDa; 5kD
Cell Pathway	Cytoplasm, cytoskeleton.
Tissue Specificity	Expressed in several hemopoietic cell lines and lymphoid malignant cells. Decreased levels in myeloma cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

