



SPAT4 rabbit pAb

Catalog No	YP-Ab-09329
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SPATA4 TSARG2
Immunogen	Synthesized peptide derived from human SPAT4 AA range: 62-112
Specificity	This antibody detects endogenous levels of human SPAT4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPAT4;SPATA4 TSARG2;Spermatogenesis-associated protein 4 (Testis and spermatogenesis cell-related protein 2) (Testis spermatocyte apoptosis-related gene 2 protein)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Highly expressed in testis, the expression is observed precisely in seminiferous tubules.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

