



AGAP5 rabbit pAb

Catalog No	YP-Ab-09332
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	AGAP5 CTGLF2
Immunogen	Synthesized peptide derived from human AGAP5 AA range: 417-467
Specificity	This antibody detects endogenous levels of human AGAP5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AGAP5;AGAP5 CTGLF2;Arf-GAP with GTPase;ANK repeat and PH domain-containing protein 5 (AGAP-5) (Centaurin-gamma-like family member 2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

