



YN4493

FCGRC rabbit pAb

Catalog No	YP-Ab-09333
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FCGR1C IGFCRC
Immunogen	Synthesized peptide derived from human FCGRC AA range: 200-250
Specificity	This antibody detects endogenous levels of human FCGRC
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FCGRC;FCGR1C IGFCRC;Putative high affinity immunoglobulin gamma Fc receptor IC (IgG Fc receptor IC) (Fc-gamma RIC) (FcRIC) (hFcgammaRIC)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images