



NECP1 rabbit pAb

Catalog No	YP-Ab-09335
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NECAP1
Immunogen	Synthesized peptide derived from human NECP1 AA range: 109-159
Specificity	This antibody detects endogenous levels of human NECP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NECP1;NECAP1;Adaptin ear-binding coat-associated protein 1 (NECAP endocytosis-associated protein 1) (NECAP-1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasmic vesicle, clathrin-coated vesicle membrane . Cell membrane . Colocalizes with AP-2 at the plasma membrane. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

