



COMD5 rabbit pAb

Catalog No	YP-Ab-09336
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	COMMD5 HT002
Immunogen	Synthesized peptide derived from human COMD5 AA range: 47-97
Specificity	This antibody detects endogenous levels of human COMD5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	COMD5;COMMD5 HT002;COMM domain-containing protein 5 (Hypertension-related calcium-regulated gene protein) (HCaRG)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Highly expressed in heart, stomach, jejunum, kidney, liver, and adrenal gland. Expression was generally higher in adult organs than in fetal tissues, particularly in heart, kidney, and liver.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

