



PHRF1 rabbit pAb

Catalog No	YP-Ab-09339
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PHRF1 KIAA1542
Immunogen	Synthesized peptide derived from human PHRF1 AA range: 1327-1377
Specificity	This antibody detects endogenous levels of human PHRF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AA673488;CTD binding SR like protein rA9;KIAA1542;PHD and RING finger domain containing protein 1;PHD and ring finger domains 1;PHRF 1;RNF221; PHRF1;PHRF1 KIAA1542;PHD and RING finger domain-containing protein 1
Observed Band	181kD
Calculated Molecular Weight	179kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

