



KLF9 rabbit pAb

Catalog No	YP-Ab-09346
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KLF9 BTEB BTEB1
Immunogen	Synthesized peptide derived from human KLF9 AA range: 49-99
Specificity	This antibody detects endogenous levels of human KLF9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KLF9;KLF9 BTEB BTEB1;Krueppel-like factor 9 (Basic transcription element-binding protein 1) (BTE-binding protein 1) (GC-box-binding protein 1) (Transcription factor BTEB1)
Observed Band	27kD
Calculated Molecular Weight	27kD
Cell Pathway	Nucleus .
Tissue Specificity	Epidermis (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

