



UNC79 rabbit pAb

Catalog No	YP-Ab-09348
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	UNC79 KIAA1409
Immunogen	Synthesized peptide derived from human UNC79 AA range: 2307-2357
Specificity	This antibody detects endogenous levels of human UNC79
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1409;Protein unc-79 homolog;UNC79;UNC79 KIAA1409
Observed Band	295 kDa
Calculated Molecular Weight	295 kDa
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

