



SPAT1 rabbit pAb

Catalog No	YP-Ab-09349
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SPATA1
Immunogen	Synthesized peptide derived from human SPAT1 AA range: 369-419
Specificity	This antibody detects endogenous levels of human SPAT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AGT;AGT1;Agxt;AGXT1;Alanine glyoxylate aminotransferase;Alanine glyoxylate aminotransferase3;Alanine--glyoxylate aminotransferase;EC 2.6.1.44;EC 2.6.1.51;Hepatic peroxisomal alanine glyoxylate aminotransferase;Hepatic peroxisomal alanine:glyoxylate aminotransferase;L alanine glyoxylate aminotransferase 1;MS773;PH1;Serine pyruvate aminotransferase;Serine--pyruvate aminotransferase;mitochondrial;Serine:pyruvate aminotransferase; SPAT;SPT;SPYA_HUMAN;TLH6;SPAT1;SPATA1;Spermatogenesis-associated protein 1 (Sperm-specific protein SP-2)
Observed Band	
Calculated Molecular Weight	43kD
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, acrosome .
Tissue Specificity	
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

