



VIR rabbit pAb

Catalog No	YP-Ab-09351
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KIAA1429 MSTP054
Immunogen	Synthesized peptide derived from human VIR AA range: 1285-1335
Specificity	This antibody detects endogenous levels of human VIR
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	VIR;KIAA1429 MSTP054;Protein virilizer homolog
Observed Band	
Calculated Molecular Weight	23 kDa
Cell Pathway	Nucleus speckle . Nucleus, nucleoplasm . Cytoplasm . Mainly nuclear with some fraction located in the cytoplasm. ZC3H13 is required to anchor component of the MACOM subcomplex, such as VIRMA, in the nucleus. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

