



SOSB1 rabbit pAb

Catalog No	YP-Ab-09352
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NABP2 OBFC2B SSB1 LP3587
Immunogen	Synthesized peptide derived from human SOSB1 AA range: 70-120
Specificity	This antibody detects endogenous levels of human SOSB1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SOSB1;NABP2 OBFC2B SSB1 LP3587;SOSS complex subunit B1 (Nucleic acid-binding protein 2) (Oligonucleotide/oligosaccharide-binding fold-containing protein 2B) (Sensor of single-strand DNA complex subunit B1) (Sensor of ssDNA subunit
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Localizes to nuclear foci following DNA damage. Foci formation is not cell-cycle dependent. Partial colocalization with RAD51 after ionizing radiation treatment.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images