



RET1G rabbit pAb

Catalog No	YP-Ab-09353
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	RAET1G
Immunogen	Synthesized peptide derived from human RET1G AA range: 192-242
Specificity	This antibody detects endogenous levels of human RET1G
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RET1G;RAET1G;Retinoic acid early transcript 1G protein
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Isoform 1]: Cell membrane ; Lipid-anchor, GPI-anchor . Endoplasmic reticulum . Mainly found intracellularly. .; [Isoform 2]: Secreted .
Tissue Specificity	Isoform 1 is highly expressed in colon and in a number of tumor cell lines and highly restricted in normal tissues. Both isoforms are frequently expressed in cell lines derived from epithelial cancers, and in primary breast cancers.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

