



F110C rabbit pAb

Catalog No	YP-Ab-09354
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM110C
Immunogen	Synthesized peptide derived from human F110C AA range: 109-159
Specificity	This antibody detects endogenous levels of human F110C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F110C;FAM110C;Protein FAM110C
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, spindle pole . Nucleus . Colocalizes with microtubules during interphase (PubMed:17499476). Detected at the mitotic spindle poles (PubMed:17499476). Colocalizes with AKT1 at the cell cortex (PubMed:19698782). .
Tissue Specificity	Detected in stomach, thyroid, trachea, adrenal gland and testis, and at low levels in prostate, ovary, intestine, colon, spinal cord and lymph node.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images