



CBLC rabbit pAb

Catalog No	YP-Ab-09358
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CBLC CBL3 RNF57
Immunogen	Synthesized peptide derived from human CBLC AA range: 219-269
Specificity	This antibody detects endogenous levels of human CBLC
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	cblC;MMACHC;RP11 291L19.3;CBLC CBL3 RNF57;Signal transduction protein CBL-C (RING finger protein 57) (SH3-binding protein CBL-3) (SH3-binding protein CBL-C)
Observed Band	32 kDa
Calculated Molecular Weight	282 aa, 32 kDa
Cell Pathway	
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

