



TEX35 rabbit pAb

Catalog No	YP-Ab-09359
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	Tex35 C1orf49
Immunogen	Synthesized peptide derived from human TEX35 AA range: 96-146
Specificity	This antibody detects endogenous levels of human TEX35
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TEX35;Tex35 C1orf49;Testis-expressed sequence 35 protein
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Testis-specific.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

