



SIRPD rabbit pAb

Catalog No	YP-Ab-09361
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SIRPD PTPNS1L2
Immunogen	Synthesized peptide derived from human SIRPD AA range: 105-155
Specificity	This antibody detects endogenous levels of human SIRPD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PTPNS1L2;SIRP delta;SIRPD;SIRPD PTPNS1L2;Signal-regulatory protein delta (SIRP-delta) (Protein tyrosine phosphatase non-receptor type substrate 1-like 2)
Observed Band	19 kDa
Calculated Molecular Weight	197 aa, 22 kDa
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

