



LHPL1 rabbit pAb

Catalog No	YP-Ab-09362
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LHFPL1 UNQ5824/PRO19643
Immunogen	Synthesized peptide derived from human LHPL1 AA range: 124-174
Specificity	This antibody detects endogenous levels of human LHPL1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LHPL1;LHFPL1 UNQ5824/PRO19643;Lipoma HMGIC fusion partner-like 1 protein
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed. Expressed at high levels in lung, thymus, skeletal muscle, colon and ovary.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

