



DC8L2 rabbit pAb

Catalog No	YP-Ab-09363
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	DCAF8L2 WDR42C
Immunogen	Synthesized peptide derived from human DC8L2 AA range: 458-508
Specificity	This antibody detects endogenous levels of human DC8L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DC8L2;DCAF8L2 WDR42C;DDB1- and CUL4-associated factor 8-like protein 2 (WD repeat-containing protein 42C)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

