



SRSF2 rabbit pAb

Catalog No	YP-Ab-09364
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SRSF2 SFRS2
Immunogen	Synthesized peptide derived from human SRSF2 AA range: 15-65
Specificity	This antibody detects endogenous levels of human SRSF2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PR264;Protein PR264;SC 35;SC35;SC-35;SRSF2;SRSF2 SFRS2;Serine/arginine-rich splicing factor 2 (Protein PR264) (Splicing component;35 kDa) (Splicing factor SC35) (SC-35) (Splicing factor;arginine/serine-rich 2)
Observed Band	25-35 kDa
Calculated Molecular Weight	25 kDa
Cell Pathway	Nucleus . Nucleus, nucleoplasm . Nucleus speckle . Phosphorylation by SRPK2 provokes its redistribution from the nuclear speckle to nucleoplasm.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

