



DPF3 rabbit pAb

Catalog No	YP-Ab-09365
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	DPF3 BAF45C CERD4
Immunogen	Synthesized peptide derived from human DPF3 AA range: 45-95
Specificity	This antibody detects endogenous levels of human DPF3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DPF3;DPF3 BAF45C CERD4;Zinc finger protein DPF3 (BRG1-associated factor 45C) (BAF45C) (Zinc finger protein cer-d4)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

