



ASCL4 rabbit pAb

Catalog No	YP-Ab-09367
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	ASCL4 BHLHA44 HASH4
Immunogen	Synthesized peptide derived from human ASCL4 AA range: 81-131
Specificity	This antibody detects endogenous levels of human ASCL4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Achaete scute complex homolog 4 (Drosophila);Achaete scute complex like 4; Achaete scute like protein 4;bHLHa44;Class II bHLH protein ASCL4;HASH4; ASCL4;ASCL4 BHLHA44 HASH4;Achaete-scute homolog 4 (ASH-4) (hASH4) (Achaete-scute-like protein 4) (Class A basic helix-loop-helix protein 44) (bHLHa44)
Observed Band	19kD
Calculated Molecular Weight	19kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in skin. 7-fold higher expression in fetal skin than in adult skin. Weak expression also detected in fetal lung, aorta and brain, and in adult stomach, kidney, ovary and breast.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

