



LRIT3 rabbit pAb

Catalog No	YP-Ab-09368
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	LRIT3
Immunogen	Synthesized peptide derived from human LRIT3 AA range: 156-206
Specificity	This antibody detects endogenous levels of human LRIT3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRIT3;Leucine-rich repeat;immunoglobulin-like domain and transmembrane domain-containing protein 3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell projection, dendrite . Perikaryon . Endoplasmic reticulum membrane ; Single-pass type I membrane protein . Punctate expression at dendrite tips. .
Tissue Specificity	Detected in the outer plexiform layer (OPL) of the retina where it localizes to ON-bipolar cells (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

