



D108A rabbit pAb

Catalog No	YP-Ab-09369
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	DEFB108P1 DEFB108A; DEFB108P2
Immunogen	Synthesized peptide derived from human D108A AA range: 22-72
Specificity	This antibody detects endogenous levels of human D108A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	D108A;DEFB108P1 DEFB108A;DEFB108P2;Putative beta-defensin 108A (Defensin;beta 108A) (Putative beta-defensin 108B pseudogene 1/2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

