



HERC1 rabbit pAb

Catalog No	YP-Ab-09371
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	HERC1
Immunogen	Synthesized peptide derived from human HERC1 AA range: 509-559
Specificity	This antibody detects endogenous levels of human HERC1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	p532;p619;HERC1;Probable E3 ubiquitin-protein ligase HERC1 (EC 6.3.2.-) (HECT domain and RCC1-like domain-containing protein 1) (p532) (p619)
Observed Band	532-550 kDa
Calculated Molecular Weight	532 kDa
Cell Pathway	Membrane; Peripheral membrane protein. Cytoplasm, cytosol. Golgi apparatus. Recruited onto actin-rich surface protrusions.
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

