



CDAC1 rabbit pAb

Catalog No	YP-Ab-09372
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CDADC1
Immunogen	Synthesized peptide derived from human CDAC1 AA range: 288-338
Specificity	This antibody detects endogenous levels of human CDAC1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CDAC1;CDADC1;Cytidine and dCMP deaminase domain-containing protein 1 (Testis development protein NYD-SP15)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Widely expressed. Expressed at high levels in the testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

