



GATS rabbit pAb

Catalog No	YP-Ab-09373
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GATS STAG3OS
Immunogen	Synthesized peptide derived from human GATS AA range: 110-160
Specificity	This antibody detects endogenous levels of human GATS
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GATS like protein 2;GATSL2;GATS;GATS STAG3OS;Putative protein GATS (STAG3 opposite strand transcript protein)
Observed Band	36 kDa
Calculated Molecular Weight	329 aa, 36 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

