



CHIC1 rabbit pAb

Catalog No	YP-Ab-09375
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CHIC1 BRX
Immunogen	Synthesized peptide derived from human CHIC1 AA range: 51-101
Specificity	This antibody detects endogenous levels of human CHIC1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CHIC1;CHIC1 BRX;Cysteine-rich hydrophobic domain 1 protein (Brain X-linked protein)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane . Cytoplasmic vesicle . Also present at a Golgi-like vesicular compartment and at scattered vesicles. .
Tissue Specificity	Equally expressed in various parts of the brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

