



YN4536

HMHB1 rabbit pAb

Catalog No	YP-Ab-09376
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	HMHB1
Immunogen	Synthesized peptide derived from human HMHB1 AA range: 1-24
Specificity	This antibody detects endogenous levels of human HMHB1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HMHB1;Minor histocompatibility protein HB-1 [Cleaved into: Minor histocompatibility antigen HB-1 (mHag HB-1)]
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in acute lymphoblastic leukemia B-cells and Epstein-Barr virus-transformed B-cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images