



WDR81 rabbit pAb

Catalog No	YP-Ab-09377
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	WDR81
Immunogen	Synthesized peptide derived from human WDR81 AA range: 525-575
Specificity	This antibody detects endogenous levels of human WDR81
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	WD repeat-containing protein 81;WDR81
Observed Band	74-80 kDa, 250 kDa
Calculated Molecular Weight	1941 aa, 212 kDa
Cell Pathway	Early endosome membrane ; Peripheral membrane protein . Late endosome membrane . Lysosome membrane . Cytoplasmic vesicle, autophagosome membrane . Mitochondrion . Cytoplasm, cytosol .
Tissue Specificity	Widely expressed. In the brain, highest levels in cerebellum and corpus callosum.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

