



PP6R1 rabbit pAb

Catalog No	YP-Ab-09380
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PPP6R1 KIAA1115 PP6R1 SAPS1
Immunogen	Synthesized peptide derived from human PP6R1 AA range: 136-186
Specificity	This antibody detects endogenous levels of human PP6R1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1115;PP6R1;PPP6R1;SAP190;SAPS domain family member 1;SAPS domain family;member 1;SAPS1;PPP6R1 KIAA1115 PP6R1 SAPS1;Serine/threonine-protein phosphatase 6 regulatory subunit 1 (SAPS domain family member 1)
Observed Band	120 kDa
Calculated Molecular Weight	881aa,97 kDa; 943aa,103 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Ubiquitous with higher expression in testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

