



KHNYN rabbit pAb

Catalog No	YP-Ab-09381
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KHNYN KIAA0323
Immunogen	Synthesized peptide derived from human KHNYN AA range: 176-226
Specificity	This antibody detects endogenous levels of human KHNYN
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein KHNYN;KHNYN;KH and NYN domain-containing protein;KHNYN KIAA0323;Protein KHNYN (KH and NYN domain-containing protein)
Observed Band	
Calculated Molecular Weight	75 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

