



BPIB6 rabbit pAb

Catalog No	YP-Ab-09383
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	BPIFB6 BPIL3
Immunogen	Synthesized peptide derived from human BPIB6 AA range: 268-318
Specificity	This antibody detects endogenous levels of human BPIB6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BPIB6;BPIFB6 BPIL3;BPI fold-containing family B member 6 (Bactericidal/permeability-increasing protein-like 3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Detected at very low levels in normal tonsils, and at higher levels in hypertrophic tonsils.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

