



GBRP rabbit pAb

Catalog No	YP-Ab-09384
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GABRP
Immunogen	Synthesized peptide derived from human GBRP AA range: 181-231
Specificity	This antibody detects endogenous levels of human GBRP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GABA(A) receptor subunit pi;GABRP;Gamma aminobutyric acid (GABA) A receptor pi;Gamma aminobutyric acid A receptor pi;Gamma aminobutyric acid receptor pi subunit;Gamma aminobutyric acid receptor subunit pi;Gamma-aminobutyric acid receptor subunit pi;GBRP_HUMAN;MGC126386;MGC126387;GBRP;Gamma-aminobutyric acid receptor subunit pi (GABA(A) receptor subunit pi)
Observed Band	
Calculated Molecular Weight	5kDa; 51kD
Cell Pathway	Cell junction, synapse, postsynaptic cell membrane; Multi-pass membrane protein. Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Most abundant in the uterus, also expressed in lung, thymus and prostate.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

