



PPM1A rabbit pAb

Catalog No	YP-Ab-09385
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PPM1A PPPM1A
Immunogen	Synthesized peptide derived from human PPM1A AA range: 284-334
Specificity	This antibody detects endogenous levels of human PPM1A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PP2C ALPHA;PP2CA;PPM1A;PPPM1A;Protein phosphatase 1A;Protein phosphatase 1A;PPM1A PPPM1A;Protein phosphatase 1A (EC 3.1.3.16) (Protein phosphatase 2C isoform alpha) (PP2C-alpha) (Protein phosphatase 1A)
Observed Band	42 kDa
Calculated Molecular Weight	382 aa, 42 kDa
Cell Pathway	Nucleus . Cytoplasm, cytosol . Membrane ; Lipid-anchor . Weakly associates at the membrane and N-myristoylation mediates the membrane localization. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

