



ENDOU rabbit pAb

Catalog No	YP-Ab-09386
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ENDO U
Immunogen	Synthesized peptide derived from human ENDOU AA range: 32-82
Specificity	This antibody detects endogenous levels of human ENDOU
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	26 serine protease;ENDO U;P11;Placental protein 11;PP11;Protein endoU; PRSS26;Poly(U)-specific endoribonuclease (EC 3.1.-.-) (Placental protein 11) (PP11) (Protein endoU) (Uridylate-specific endoribonuclease)
Observed Band	35-40 kDa
Calculated Molecular Weight	369 aa, 42 kDa
Cell Pathway	Secreted .
Tissue Specificity	Placental-specific, but also associated with various malignant neoplasms.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

