



MTG8R rabbit pAb

Catalog No	YP-Ab-09387
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CBFA2T2 EHT MTGR1
Immunogen	Synthesized peptide derived from human MTG8R AA range: 67-117
Specificity	This antibody detects endogenous levels of human MTG8R
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AML1 MTG8;cbfa2t2;Core binding factor runt domain alpha subunit 2 translocated to 2;core-binding factor;runt domain;alpha subunit 2; DKFZp313F2116;EHT;ETO homolog on chromosome 2;ETO homologous on chromosome 2;MTG 8 like protein;MTG 8 related protein 1;MTG8 like protein; MTG8 related protein 1;MTG8-like protein;MTG8-related protein 1; MTG8R_HUMAN;MTGR 1;Myeloid translocation gene related protein 1;Myeloid translocation related protein 1;Myeloid translocation-related protein 1;p85; Protein CBFA2T2;ZMYND 3;ZMYND3;MTG8R;CBFA2T2 EHT MTGR1;Protein CBFA2T2 (ETO homologous on chromosome 20) (MTG8-like protein) (MTG8-related protein 1) (Myeloid translocation-related protein 1) (p85)
Observed Band	
Calculated Molecular Weight	8 kDa, 65 kDa; 67kD
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitously expressed in fetal and adult tissues. Highly expressed in adult brain, heart, lung, kidney, lymph node, appendix, thymus, testis, uterus, small intestine, prostate and thymus.
Function	



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

