



S35B3 rabbit pAb

Catalog No	YP-Ab-09388
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SLC35B3 C6orf196 PAPST2 CGI-19
Immunogen	Synthesized peptide derived from human S35B3 AA range: 258-308
Specificity	This antibody detects endogenous levels of human S35B3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	S35B3;SLC35B3 C6orf196 PAPST2 CGI-19;Adenosine 3'-phospho 5'-phosphosulfate transporter 2 (3'-phosphoadenosine 5'-phosphosulfate transporter) (PAPS transporter 2) (Solute carrier family 35 member B3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein .
Tissue Specificity	Preferentially and highly expressed in colon.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

