



MREG rabbit pAb

Catalog No	YP-Ab-09390
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MREG DSU HDCGA21P
Immunogen	Synthesized peptide derived from human MREG AA range: 98-148
Specificity	This antibody detects endogenous levels of human MREG
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DSU;melanoregulin;MREG;WDT2;MREG DSU HDCGA21P;Melanoregulin (Dilute suppressor protein homolog)
Observed Band	26 kDa
Calculated Molecular Weight	25kDa
Cell Pathway	Apical cell membrane ; Peripheral membrane protein . Melanosome membrane ; Lipid-anchor . Lysosome membrane ; Lipid-anchor . Cytoplasmic vesicle membrane . Localizes to the inner segment and basal outer segment of rods in the retina. .
Tissue Specificity	Expressed in photoreceptor cells (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images