



VP33A rabbit pAb

Catalog No	YP-Ab-09392
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VPS33A
Immunogen	Synthesized peptide derived from human VP33A AA range: 272-322
Specificity	This antibody detects endogenous levels of human VP33A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	VP33A;VPS33A;Vacuolar protein sorting-associated protein 33A (hVPS33A)
Calculated Molecular Weight	66kD
Cell Pathway	Cytoplasmic vesicle . Late endosome membrane ; Peripheral membrane protein ; Cytoplasmic side . Lysosome membrane ; Peripheral membrane protein ; Cytoplasmic side . Early endosome . Cytoplasmic vesicle, autophagosome . Cytoplasmic vesicle, clathrin-coated vesicle .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images