



SRB4D rabbit pAb

Catalog No	YP-Ab-09393
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SRCRB4D
Immunogen	Synthesized peptide derived from human SRB4D AA range: 252-302
Specificity	This antibody detects endogenous levels of human SRB4D
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SRB4D;SRCRB4D;Scavenger receptor cysteine-rich domain-containing group B protein (S4-SRCRB)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

