



AFF2 rabbit pAb

Catalog No	YP-Ab-09395
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	AFF2 FMR2 OX19
Immunogen	Synthesized peptide derived from human AFF2 AA range: 1104-1154
Specificity	This antibody detects endogenous levels of human AFF2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AF4/FMR2 family member 2;FMR2;FMR2P;FRAXE;MRX2;AFF2;AFF2 FMR2 OX19;AF4/FMR2 family member 2 (Fragile X E mental retardation syndrome protein) (Fragile X mental retardation 2 protein) (FMR2P) (Protein FMR-2) (Protein Ox19)
Observed Band	145 kDa
Calculated Molecular Weight	1311 aa, 145 kDa
Cell Pathway	Nucleus speckle . When splicing is inhibited, accumulates in enlarged speckles.
Tissue Specificity	Brain (most abundant in hippocampus and amygdala), placenta and lung.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

