



DPP3 rabbit pAb

Catalog No	YP-Ab-09396
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DPP3
Immunogen	Synthesized peptide derived from human DPP3 AA range: 152-202
Specificity	This antibody detects endogenous levels of human DPP3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Dipeptidyl aminopeptidase III;Dipeptidyl arylamidase III;dipeptidyl peptidase 3; Dipeptidyl peptidase III;DPP III;DPP3;DPPIII;Dipeptidyl peptidase 3 (EC 3.4.14.4) (Dipeptidyl aminopeptidase III) (Dipeptidyl arylamidase III) (Dipeptidyl peptidase III) (DPP III)
Observed Band	69 kDa, 83 kDa
Calculated Molecular Weight	83 kDa
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Detected in placenta (at protein level) (PubMed:3233187). Detected in erythrocytes (at protein level) (PubMed:1515063).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

