



ENTP2 rabbit pAb

Catalog No	YP-Ab-09398
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ENTPD2 CD39L1
Immunogen	Synthesized peptide derived from human ENTP2 AA range: 390-440
Specificity	This antibody detects endogenous levels of human ENTP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NTPDase 2;NTPDase2;CD39 antigen like 1;CD39 antigen-like 1;CD39 like1; CD39L1;CD39like1;ecto ATP diphosphohydrolase 2;ecto ATPase 2;ecto ATPDase 2;Ecto-ATP diphosphohydrolase 2;Ecto-ATPase 2;Ecto-ATPDase 2; ectoATPase 2;ectoATPDase 2;Ectonucleoside triphosphate diphosphohydrolase 2;ENTP2_HUMAN;Entpd2;ENTP2;ENTPD2 CD39L1;Ectonucleoside triphosphate diphosphohydrolase 2 (NTPDase 2) (EC 3.6.1.-) (CD39 antigen-like 1) (Ecto-ATP diphosphohydrolase 2) (Ecto-ATPDase 2) (Ecto-ATPase 2)
Observed Band	
Calculated Molecular Weight	54 kDa; 54kD(
Cell Pathway	[Isoform Long]: Cell membrane; Multi-pass membrane protein.; [Isoform Short]: Endoplasmic reticulum membrane; Multi-pass membrane protein.; [Isoform gamma]: Endoplasmic reticulum membrane; Multi-pass membrane protein.
Tissue Specificity	Brain, placenta, skeletal muscle, kidney, pancreas, heart, ovary, testis, colon, small intestine, prostate and pancreas. No expression in adult thymus, spleen, lung, liver and peripheral blood leukocytes.
Function	
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images