



CNOTA rabbit pAb

Catalog No	YP-Ab-09399
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CNOT10
Immunogen	Synthesized peptide derived from human CNOTA AA range: 20-70
Specificity	This antibody detects endogenous levels of human CNOTA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCR4 NOT transcription complex subunit 1;CCR4-NOT transcription complex subunit 1;CNOT 1;CNOT1;CNOTA_HUMAN;DKFZp434K115;FLJ1289;FLJ13165;CNOTA;CNOT10;CCR4-NOT transcription complex subunit 10
Observed Band	
Calculated Molecular Weight	82 kDa; 82kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

