



ZN444 rabbit pAb

Catalog No	YP-Ab-09400
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF444 EZF2 ZSCAN17
Immunogen	Synthesized peptide derived from human ZN444 AA range: 189-239
Specificity	This antibody detects endogenous levels of human ZN444
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Endothelial zinc finger protein 2;EZF 2;EZF-2;EZF2;FLJ111137;ZCCAN17;Zinc finger and SCAN domain-containing protein 17;Zinc finger protein 444; ZN444_HUMAN;ZNF444;ZSCAN17;ZN444;ZNF444 EZF2 ZSCAN17;Zinc finger protein 444 (Endothelial zinc finger protein 2) (EZF-2) (Zinc finger and SCAN domain-containing protein 17)
Observed Band	
Calculated Molecular Weight	35kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

